CE  
0120


# BIB<sup>®</sup>

## (Balloon in Balloon)

### DILATATION CATHETER

Instructions for Use.....	2
Mode d'emploi .....	7
Istruzioni per l'uso.....	12
Gebrauchsanweisung.....	17
Instrucciones de utilización.....	22
Bruksanvisning .....	27
Brugsanvisning .....	32
Gebruiksaanwijzing .....	37
Instruções de utilização .....	42
Kullanım Talimatı .....	47

**CAUTION:** Federal (USA) Law restricts this device to sale by or on the order of a physician.

## Instructions for Use:

### INDICATIONS

Recommended for Percutaneous Transluminal Angioplasty (PTA) of the femoral, iliac and renal arteries. These catheters are not designed to be used in the coronary arteries.

**The safety and effectiveness of this device for use as a stent delivery system in the placement of intravascular stents has not been established.**

### DESCRIPTION

The NuMED BIB<sup>®</sup> (Balloon in Balloon) Catheter is a triaxial design catheter. Two lumens are used to inflate the balloons while one lumen is for tracking over a guidewire. The radiopaque platinum marker(s) are placed beneath the 'working area' of the balloon. The inner balloon is ½ of the outer balloon diameter and 1 cm shorter. Each balloon inflates to the stated diameter and length at a specific pressure. The balloon size is ±10% at rated burst pressure (RBP). The RBP is different for each size. Check the package label for the RBP. It is important that the balloon not be inflated beyond the RBP.

The purpose of the double balloon catheter is to apply an incremental inflation for the purpose of opening an artery. The inner balloon provides initial expansion of the artery and also acts as a tool to hold the catheter in the artery while the outer balloon is inflated. The outer balloon is then inflated providing the remainder of the expansion.

### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas. Sterile and non-pyrogenic if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

### WARNING

- Do not exceed the RBP. An inflation device with pressure gauge is recommended to monitor pressure. Pressure in excess of the RBP can cause balloon rupture and potential inability to withdraw the catheter through the introducer sheath.
- Confirm that the distal end of the introducing sheath is at least 2.5cm back from the most proximal image markers before inflating the outer balloon. Failure to do so may stretch the outer tubing and severely hinder balloon deflation.
- In PTA, the dilated balloon should not markedly exceed the diameter of the artery lying just proximal to the stenosis.
- Use only appropriate balloon inflation medium. Do not use air or gaseous medium to inflate the balloon.
- Use two appropriate size inflation devices with pressure gauges for inflation.
- Do not advance the guidewire, balloon dilatation catheter, or any other component if resistance is met, without first determining the cause and taking remedial action.
- This catheter is not recommended for pressure measurement or fluid injection.
- Do not remove the guidewire from the catheter at any time during the procedure except when the procedure has been completed.
- This device is intended for single use only. Do not resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of cross-contamination.

### PRECAUTIONS

- Dilatation procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- Guidewires are delicate instruments. Care should be exercised while handling to help prevent the possibility of breakage.

- Careful attention must be paid to the maintenance of tight catheter connections by aspiration before proceeding to avoid air introduction into the system.
- Under no circumstances should any portion of the catheter system be advanced against resistance. The cause of the resistance should be identified with fluoroscopy and action taken to remedy the problem.
- If resistance is felt upon removal, then the balloon, guidewire and the sheath should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
- Before removing the catheter from the sheath, it is very important that the balloon is completely deflated.
- Proper functioning of the catheter depends on its integrity. Care should be used when handling the catheter. Damage may result from kinking, stretching, or forceful wiping of the catheter.

## POTENTIAL COMPLICATIONS

- Potential complications related to the introduction of the catheter into the body include, but are not limited to, the following: infection, air embolism, and hematoma formation.
- Potential balloon separation following balloon rupture or abuse and the subsequent need to use a snare or other medical interventional techniques to retrieve the pieces.
- Complications associated with PTA include, but are not limited to: clot formation and embolism, nerve damage, vascular perforation requiring surgical repair, damage to the vascular intima, cerebral accident, cardiac arrhythmias, myocardial infarction, or death. For specifics, refer to: Fellows, K. et al.: Acute Complications of Catheter therapy for Congenital Heart Disease, Amer Journ of Cardiol, 60;679(1987).

**NOTE:** There have been infrequent reports of larger diameter balloons bursting circumferentially, possibly due to a combination of tight focal strictures in large vessels. In any instance of a balloon rupture while in use, it is recommended that a sheath be placed over the ruptured balloon prior to withdrawal through the entry site. This can be accomplished by cutting off the proximal end of the catheter and slipping an appropriately sized sheath over the catheter into the entry site. For specific technique, refer to: Tegtmeier, Charles J., M.D. & Bezirdijan Diran R., M.D. "Removing the Stuck, Ruptured Angioplasty Balloon Catheter." Radiology, Volume 139, 231-232, April 1981.

## INSPECTION AND PREPARATION

1. Using proper sterile technique, open the catheter package and remove BIB<sup>®</sup> catheter. Inspect the catheter for kinks prior to use. **DO NOT REMOVE BALLOON PROTECTOR FROM BALLOON.** The catheter lumen (green hub), the inner balloon (indigo hub), and the outer balloon (orange hub) are appropriately labeled. Flush catheter lumen (green hub) with heparinized flush and insert .035 guidewire. **DO NOT ATTEMPT TO PURGE BALLOONS WITHOUT A GUIDEWIRE THROUGH THE CATHETER LUMEN.**
2. Prepare an inflation solution of 40% by volume of contrast medium (such as Renografin<sup>®</sup> 76%, Omnipaque<sup>®</sup> 300 etc.) and flush solution.
3. Fill and purge air from the inflation devices (2 needed for BIB<sup>®</sup> catheter). Rotate stopcock handle to "off" position relative to inflation device. Apply negative pressure with the inflation device and tip to remove air from manometer. Repeat as necessary. Attach one inflation device to each inflation port of the BIB<sup>®</sup> catheter using a three way stopcock with rotating adapter. The indigo hub is the smaller inner balloon, the orange hub is the larger outer balloon.
4. Purging of the two balloons is done with negative pressure only. **DO NOT INFLATE EITHER BALLOON PRIOR TO USE.**
5. Attach a syringe capable of maintaining negative pressure (e.g. VacLok Syringe<sup>®</sup> from Merit Medical or Becton-Dickinson<sup>®</sup> luer lock Syringe) to the other port of the stopcock. With the stopcock handle in the "off" position to the syringe, apply negative pressure to the syringe and lock. **ALWAYS START WITH THE INNER BALLOON (INDIGO HUB).** Turn the handle of the stopcock to apply negative pressure to the balloon. Tap the catheter shaft to facilitate movement of bubbles, and then shut off the negative pressure. Some contrast will passively move from inflation device into the balloon. Repeat, alternating negative pressure and

- passive contrast, until you no longer have bubbles when applying negative pressure. It will be necessary to repeat this process several times. Repeat process with the outer balloon.
6. Leave both balloons on continuous negative pressure for insertion. Remove balloon protector prior to removing the guidewire. Flush catheter lumen with heparinized flush.

## **INSTRUCTIONS FOR USE**

Prior to angioplasty, carefully examine all equipment to be used during the procedure, including the catheter, to verify proper function and that the catheter size is suitable for the specific procedure for which it is intended. Also, inflate the dilatation catheter to the appropriate RBP and deflate to verify proper function.

1. Enter the vessel percutaneously using the standard Seldinger technique over the appropriate guidewire for the size catheter being used.
2. Under fluoroscopic guidance advance the guidewire to the desired position. Pass the catheter over the guidewire. An introducer should be utilized to facilitate catheter insertion.
3. Advance the catheter across the lesion with fluoroscopic guidance using accepted percutaneous transluminal angioplasty technique (see references). In most patients, the balloon should meet with minimal resistance to insertion. Do not advance the catheter unless the guidewire is in place. (Place the catheter to position the mid-length of the balloon within the lesion. Radiopaque bands define both the inner and outer balloon shoulders.)
4. The distal lumen is provided for guidewire tracking. An inflation device with pressure gauge is required to monitor inflation pressure [refer to package label for RBP].
5. Inflate the inner balloon to RBP or lower and close stopcock. Immediately inflate the outer balloon to RBP to achieve dilatation. DO NOT exceed the RBP.
6. Immediately after dilatation, deflate both balloons. Both balloons are deflated simultaneously. Deflate the balloon by drawing a vacuum with an inflation device with pressure gauge. Note: The greater the vacuum applied and held during withdrawal, the lower the deflated balloon profile. Gently withdraw the catheter. As the balloon exits the vessel, use a smooth, gentle, steady motion. If resistance is felt upon removal, then the balloon, guidewire and the sheath should be removed together as a unit under fluoroscopic guidance, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
7. Apply pressure to the insertion site according to standard practice or hospital protocol for percutaneous vascular procedures.

## **WARNING:**

NuMED catheters are placed in the extremely hostile environment of the human body. Catheters may fail to function for a variety of causes including, but not limited to, medical complications or failure of catheters by breakage. In addition, despite the exercise of all due care in design, component selection, manufacture and testing prior to sale, catheters may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of catheters will not occur or that the body will not react adversely to the placement of catheters or that medical complications will not follow the use of catheters. NuMED cannot warrant or guarantee NuMED accessories because the structure of the accessories may be damaged by improper handling before or during use. Therefore, no representations or warranties are made concerning them.

## **Warranty and Limitations**

Catheters and accessories are sold in an 'as is' condition. The entire risk as to the quality and performance of the catheter is with the buyer. NuMED disclaims all warranties, expressed or implied, with respect to catheters and accessories, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. NuMED shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind NuMED to any representation or warranty with respect to catheters and accessories.

## BIB<sup>®</sup> Balloon Sizing Chart

Applied Press.	Inner Balloon 4.0 (mm)	Outer Balloon 8.0 (mm)	Inner Balloon 5.0 (mm)	Outer Balloon 10.0 (mm)	Inner Balloon 6.0 (mm)	Outer Balloon 12.0 (mm)	Inner Balloon 7.0 (mm)	Outer Balloon 14.0 (mm)	Inner Balloon 7.0 (mm)	Outer Balloon 15.0 (mm)
1.0 ATM	3.24	6.67	3.72	8.05	5.11	10.17	6.17	12.14	6.17	13.25
2.0 ATM	3.41	6.84	4.11	8.40	5.32	10.54	6.36	12.49	6.41	13.71
3.0 ATM	3.58	6.99	4.43	8.70	5.51	11.01	6.54	13.00	6.59	14.33
4.0 ATM	3.78	7.19	4.66	9.07	5.74	11.51	6.71	13.51	6.79	14.75
4.5 ATM										
5.0 ATM	<b>3.98</b>	7.37	<b>4.95</b>	9.35	<b>5.97</b>	11.81	<b>6.97</b>	13.88	<b>6.96</b>	<b>15.15</b>
6.0 ATM		7.56		9.61		12.06		<b>14.15</b>		
7.0 ATM		7.71		9.80		<b>12.24</b>				
8.0 ATM		7.84		9.93						
9.0 ATM		7.95		<b>10.09</b>						
10.0 ATM		<b>8.06</b>								

Applied Press.	Inner Balloon 8.0 (mm)	Outer Balloon 16.0 (mm)	Inner Balloon 9.0 (mm)	Outer Balloon 18.0 (mm)	Inner Balloon 10.0 (mm)	Outer Balloon 20.0 (mm)	Inner Balloon 11.0 (mm)	Outer Balloon 22.0 (mm)	Inner Balloon 12.0 (mm)	Outer Balloon 24.0 (mm)
1.0 ATM	6.16	14.17	7.78	15.33	8.88	17.66	9.77	19.55	10.08	22.09
2.0 ATM	6.84	14.87	8.08	16.59	9.17	18.63	10.09	20.70	10.50	23.09
3.0 ATM	7.24	15.60	8.41	17.45	9.43	19.64	10.36	<b>21.98</b>	11.09	<b>24.11</b>
4.0 ATM	7.71	16.21	8.69	<b>18.11</b>	9.63	<b>20.46</b>	10.78		11.62	
4.5 ATM							<b>11.00</b>		<b>11.96</b>	
5.0 ATM	<b>7.96</b>	<b>16.58</b>	<b>8.98</b>		<b>9.93</b>					
6.0 ATM										
7.0 ATM										
8.0 ATM										
9.0 ATM										
10.0 ATM										

The figures in bold face represent the balloon diameter @ Rated Burst Pressure.

**FOR ALL NUMED CATHETERS AN INFLATION DEVICE  
WITH PRESSURE GAUGE SHOULD BE USED.**

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
## REFERENCES:

- Athanasoulos, C.A.; PTA: General Principles, Amer Journ of Roent, 135; 893-900, (1980)
- Gruntzig A.; Kumpe, D.; Technique of PTA with the Gruntzig Balloon Catheter, Amer Journ of Roent 132:547-552 (1979)
- Katzen, B.T.; Chang, J.: PTA with the Gruntzig Balloon Catheter, Radiology, 130:623-626 (1979)
- Zeitler, E.; Gruntzig, A.; Schoop, W.; (Eds) Percutaneous Vascular Recanalization, Springer-Verlag, Heidelberg/New York (1978)
- Stanson A.W.; A Perspective of PTA, Cardiovascular Clinics, 12(2): 245-259 (1983)
- Wilson, A.R.; Fuchs, J.C.A.; PTA, Surgical Clinics of North America, 64,1121-1150, (1984)
- Tegtmeyer, C.J.; Kofler T.J.; Ayers C.A.; Renal Angioplasty: Current Status, AJR 142: 17-21, (1984)
- Sos, T.A. et. al.; The Current Role of Renal Angioplasty in the Treatment of Renovascular Hypertension, Urology Clinics of North America Vol. 11:3, 503-513, (1984)
- Fellows, K. et al.: Acute Complications of Catheter Therapy for Congenital Heart Disease, Amer Journ of Cardiol, 60:679 (1987)
- Rocchini, A.: Balloon Angioplasty for Peripheral Pulmonary Artery Stenosis, in Rao, P.S. (Ed) Transcatheter therapy in Pediatric Cardiology, Wiley-Liss, New York, PP 213-228 (1993)

	Do Not Reuse Ne pas réutiliser Non riutilizzare Nicht wiederverwenden No reutilizar	Får ej återvändas Må ikke genbruges Niet hergebruiken Não Reutilizar Tekrar Kullanmayın
	Catalogue Number Número de catalogue Numero di catalogo Katalognummer Número de catálogo	Katalognummer Katalognummer Catalogusnummer Número de Catálogo Katalog Numarasi
	Batch Code Code du lot Codice del lotto Chargencode Código del lote	Sändningskod Partikode Partijcode Código do Lote Parti Kodu
	Non-Sterile Non sterile Non sterile Nicht steril No estéril	Ikke steril Ikke-steril Niet-steriel Não Estéril Steril Değildir
	Attention: See Instructions For Use Attention : Consulter le mode d'emploi Attenzione – Consultare le Istruzioni per l'uso Achtung: Siehe Gebrauchsanweisung Atención: consultar instrucciones de uso	OBS! Se bruksanvisningen Bemærk: Se brugsanvisningen Let op: zie gebruiksaanwijzing Atenção: Consulte as instruções de utilização Dikkat: Kullanım Talimatlarını Bakın
	Keep Away From Sunlight Garder à l'abri de la lumière du soleil Proteggere dai raggi solari Sonneneinstrahlung vermeiden Mantener alejado de la luz solar	Exponera inte för direkt solljus Holdes væk fra sollys Uit direct zonlicht houden Manter Afastado da Luz Solar Güneş Işığından Uzak Tutun
	Temperature Limitation Limite de température Limitazione di temperatura Temperaturbegrenzung Limitación de la temperatura	Temperaturbegränsning Temperaturbegrænsning Temperaturlæperking Limite de Temperatura Isı Sınırlaması
	Manufacturer / Date of Manufacture Fabricant/Date de fabrication Produttore/Data di fabbricazione Hersteller/Herstellungsdatum Fabricante/Fecha de fabricación	Tillverkare/Tillverkningsdatum Producent/fremstillingsdato Fabrikant / Productiedatum Fabricante / Data de Fabrico Üretici / Üretim Tarihi
	Use By À utiliser avant Utilizzare entro Verwendbar bis Usar antes de	Använd senast Anvendes inden Gebruiken voor Usar Até Son Kullanma Tarihi
	Sterilized Using Ethylene Oxide Stérilisé à l'oxyde d'éthylène Sterilizzato con ossido di etilene Sterilisiert mit Ethylenoxid Esterilizado con óxido de etileno	Steriliserad med etylenoxid Steriliseret med ethylenoxid Gesteriliseerd met ethyleenoxide Esterilizado Usando Óxido de Etileno Etilen Oksit Kullanılarak Sterilize Edilmiştir
	Authorised Representative in the European Community Représentant autorisé dans la Communauté européenne Rappresentante autorizzato nella Comunità Europea Autorisierter Vertreter in der Europäischen Gemeinschaft Representante autorizado en la Comunidad Europea	Autoriserad representant inom EU Autoriseret repræsentant i Det Europæiske Fællesskab Gevolmachtigd vertegenwoordiger in de Europese Gemeenschap Representante Autorizado na Comunidade Europeia Avrupa Topluluğundaki Yetkili Temsilci

IFU-420CE-18

NuMED, Inc.  
2880 Main Street  
Hopkinton, New York USA 12965  
Telephone: (315) 328-4491  
Facsimile: (315) 328-4941  
email: [info@numedusa.com](mailto:info@numedusa.com)  
[www.numedforchildren.com](http://www.numedforchildren.com)

European Representative:  
G. van Wageningen B.V.  
Hallenweg 40, 5683 CT Best,  
The Netherlands  
TEL: +31-499-377388  
FAX: +31-499-377456  
email: [info@heartmedical.nl](mailto:info@heartmedical.nl)