
# NuDEL™

## COARCTATION OF THE AORTA

|                                   |    |
|-----------------------------------|----|
| Instructions for Use.....         | 2  |
| Mode d'emploi .....               | 8  |
| Istruzioni per l'uso.....         | 14 |
| Gebrauchsanweisung.....           | 21 |
| Instrucciones de utilización..... | 27 |
| Bruksanvisning .....              | 33 |
| Brugsanvisning .....              | 39 |
| Gebruiksaanwijzing .....          | 45 |
| Instruções de utilização .....    | 52 |
| Kullanım Talimatı .....           | 58 |

CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.

**INSTRUCTIONS FOR USE:****INDICATIONS**

Indicated for implantation in the native and/or recurrent coarctation of the aorta on patients with the following clinical conditions:

- Stenosis of the aorta resulting in significant anatomic narrowing as determined by angiography or non-invasive imaging, i.e. echocardiography, magnetic resonance imaging (MRI), CT scan;
- Stenosis of the aorta resulting in hemodynamic alterations, resulting in systolic pressure gradient, systemic hypertension or altered left ventricular function;
- Stenosis of the aorta where balloon angioplasty is ineffective or contraindicated;
- Stenosis diameter <20% of the adjacent vessel diameter. Stenosis that would present increased risk of vascular damage or disruption; or aneurysm associated with coarctation of the aorta.

**DESCRIPTION**

The NuDEL System includes a triaxial balloon in balloon designed catheter with a Covered CP Stent mounted on it; which is then covered by a sheath, as an all-in-one system. The catheter tip is tapered permitting safe entry of the access vessel.

The Covered CP Stent is balloon expandable and intended for permanent implant. The Covered CP Stent is composed of heat treated 90% platinum/10% iridium wire that is arranged in laser welded rows with a "zig" pattern. The number of zigs in a row can be varied and will impact the strength of the stent as well as the eventual expanded diameter and percent stent shortening, while the number of rows will determine the unexpanded length of the stent. The Covered CP Stent has an ePTFE covering attached to the stent framework. This covering acts as a fluid barrier creating a fluid tight conduit through the stent length.

The catheter portion of this system is a balloon in balloon triaxial catheter. Two lumens are used to inflate the balloons while one lumen is for tracking over a guidewire. The radiopaque markers are placed beneath the "working area" of the outer balloon. The inner balloon is ½ of the outer balloon diameter and 1 cm shorter. Each balloon inflates to the stated diameter and length at a specific pressure. The balloon size is ± 10% at rated burst pressure (RBP). The RBP is different for each size. Check the package label for RBP. It is important that the balloon not be inflated beyond the RBP. The inner balloon provides an initial expansion of the stent and acts as a tool to hold the stent in place, and if necessary reposition the stent, before outer balloon is inflated to deploy the stent, securing the stent against the vessel wall.

**HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas. Sterile and non-pyrogenic if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

**CONTRAINDICATIONS**

- Patients too small to allow safe delivery of the stent without compromise to the systemic artery used for delivery;
- Unfavorable aortic anatomy that does not dilate with high pressure balloon angioplasty;
- Occlusion or obstruction of systemic artery precluding delivery of the stent;
- Clinical or biological signs of infection;
- Active endocarditis;
- Known allergy to aspirin, other antiplatelet agents, or heparin;
- Pregnancy.

**WARNINGS**

- The sheath must be flushed with heparinised saline via the proximal side port prior to introducing the delivering system in the body.
- As with any type of implant, infection secondary to contamination of the stent may lead to aortitis, or abscess.

- The platinum/iridium stent may migrate from the site of implant.
- Over-stretching of the artery may result in rupture or aneurysm formation.
- The inflated diameter of the stent should at least equal the diameter of the intended implant site.
- Excessive force while crimping may weaken welds of the stent.
- Crimping the 8 zig stent on a balloon catheter smaller than 12mm, may cause damage to the stent.
- Excessive handling and manipulation of the covering while crimping the stent may cause the covering to tear off of the stent.
- Crimping the device in the opposite direction of the folds in the covering may cause the covering to catch while inserting into the hemostasis valve tool and introducer. This could cause the covering to tear off of the stent.
- Retracting the covered stent back in to the sheath may cause the covering to catch and tear off of the stent.
- Do not exceed the RBP. An inflation device with a pressure gauge is recommended to monitor pressure. Pressure in excess of the RBP can cause balloon rupture and potential inability to withdraw the catheter into the sheath.
- Confirm that the distal end of the outer sheath is at least 2.5 cm back from the proximal image marker before inflating the outer balloon. Failure to do so may stretch the outer tubing and severely hinder inner balloon deflation.
- It is recommended to use two appropriate size inflation devices with pressure gauges for inflation.
- Do not advance the guidewire, the combined balloon catheter in the sheath, or any other component if resistance is encountered, without first determining the cause and taking remedial action.
- This catheter is not recommended for pressure measurement.
- Do not remove the guidewire from the sheath-catheter at any time during the procedure except when the procedure has been completed.
- This device is intended for single use only. Do not resterilize and/or reuse it, as this can potentially compromise device integrity and performance, and increased risk of cross contamination and infection.

## **PRECAUTIONS**

- Use of an inflation device with pressure gauge is highly recommended during this procedure.
- Sheath-catheter manipulation and stent deployment must be conducted under fluoroscopic guidance with appropriate radiographic equipment.
- Stents are delicate devices. Exercise caution when handling the stent to prevent breakage.
- Guidewires should be handled with care to avoid kinking or breaking. When advancing the NuDEL system over the wire, the tip of the wire must be controlled at all times.
- The NuDEL system, especially at the stent, is rigid and may make negotiation through vessels difficult.
- Maintain tight catheter connections at all times.
- De-air the sheath with heparinized saline prior to insertion in the patient. Apply repeated negative pressure cycles to the balloons to replace the air in the balloons with fluid: this enhances the efficiency of inflation and avoids air introduction into the circulation in the unusual case of balloon rupture.
- The inflation diameter of the balloon used during stent delivery should approximate the diameter of the obstructive vessel and the intended implant site.
- Under no circumstances should any portion of the sheath-catheter system be advanced or removed against resistance. Use fluoroscopy to identify and resolve the resistance.
- If resistance is encountered upon removal, the whole system (balloon, guidewire and sheath) should be removed as a single unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping and withdrawing the sheath and catheter, using a gentle anticlockwise twisting motion with traction.
- The balloons must be completely deflated before retracting into the sheath.

- Proper functioning of the sheath-catheter depends on its integrity. Exercise caution when handling the catheter. Damage may result from kinking, stretching, or forceful wiping of the catheter.

## POTENTIAL COMPLICATIONS/ADVERSE EFFECTS

NOTE: Circumferential tear of the delivery balloon catheter prior to complete expansion of the stent may cause the balloon to become tethered to the stent, requiring surgical removal. In case of rupture of an adequately sized balloon after stent expansion, it can be withdrawn and a new balloon catheter exchanged over a guidewire to complete expansion of the stent.

Cardiac catheterization carries certain risks. In addition, potential complications, and related adverse effects associated with implants include, but are not limited to:

- Femoral artery injury, thrombosis or pseudoaneurysm
- Stent Migration
- Stent Fracture
- Aortic Rupture/Tear
- Hematoma
- Thrombosis/Thromboembolism
- Death
- Endocarditis
- Cell necrosis at the site of implant
- Stent Stenosis
- Aortic Aneurysm/Pseudoaneurysm
- Stent Malposition
- Sepsis/Infection
- AV fistula formation
- Transitory arrhythmia
- Bleeding
- Cerebrovascular Incident

## MRI SAFETY INFORMATION

Nonclinical testing and modeling has demonstrated that the CP Stent is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3 T
- Maximum spatial gradient magnetic field of 2500 gauss/cm (25 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning (Normal Operating Mode)

Based on nonclinical testing and modeling, under the scan conditions defined above, the CP Stent is expected to produce a maximum in vivo temperature rise of less than 2°C after 15 minutes of continuous scanning.

MR image quality may be compromised if the area of interest is in the same area, or relatively close to the position of the device. In nonclinical testing, the image artifact caused by the device extends approximately 3 mm from the CP Stent when imaged with a spin echo pulse sequence and 6 mm when imaged with a gradient echo pulse sequence and a 3 T MRI System. The lumen of the device was obscured.

The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters.

## INSPECTION AND PREPARATION

1. Using proper sterile technique, open the catheter package and remove the NuDEL system. Inspect the sheath and catheter for kinks prior to use. There will be four hubs on the system, and they are: catheter lumen (green hub), the inner balloon (indigo hub), the outer balloon (orange hub), and the outer sheath (white hub). Flush the catheter lumen (green hub) with heparinized saline and insert 0.035" guidewire. **DO NOT ATTEMPT TO PURGE BALLOONS WITHOUT A GUIDEWIRE THROUGH THE CATHETER LUMEN.**
2. Prepare an inflation solution of 30-50% by volume of contrast medium and flush solution.
3. Advance the catheter out of the sheath to expose the balloons without fully exposing the covered stent. Fill the inflation device with inflation solution. Attach one inflation device to each port of the NuDEL system using a three-way stopcock with rotating adapter. Purge air from the inflation devices (2 needed for the NuDEL system). Attach a syringe capable of maintaining negative pressure to the other port of the stopcock. **ALWAYS START WITH THE INNER BALLOON (INDIGO HUB).** Open the stopcock to the balloon and apply negative pressure with the inflation device to remove air from the balloon. Tap the catheter shaft and the inflation device repeatedly to facilitate movement of bubbles out

of the balloon and inflation device. Close the stopcock to the balloon, and expel the air into the syringe. Tipping the inflation device downward, release the negative pressure, and, open the stopcock to the balloon allowing the inflation solution to flow into balloon and its lumen. Repeat, alternating negative pressure and filling the balloon and its lumen with inflation solution until bubbles are no longer visible in the inflation device. It is necessary to repeat this process several times. Repeat the process with the outer balloon. Flushing the sheath (white hub) continuously with heparinised saline, the sheath is advanced to fully cover the mounted stent ensuring the delivering sheath is devoid of air.

4. Purging of the two balloons is done with negative pressure only. **DO NOT INFLATE EITHER BALLOON PRIOR TO USE.**
5. Leave both balloons on continuous negative pressure for insertion. Remove the 0.035" guidewire. Flush the catheter lumen with heparinised saline.

## **INSTRUCTIONS FOR USE**

### **Select Stent Size**

1. Measure the length of the target stricture to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the stricture.
2. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent.  
Note: Should more than one stent be required, place the stent most distal from the puncture site first, followed by placement of the proximal stent in tandem.
3. Measure the diameter of the reference stricture and vessel proximal and distal to the target lesion to determine the appropriate size stent and delivery system.
4. Consideration should be given to the stent shortening during dilation. See Stent Foreshortening chart.

### **Stent Deployment**

1. When preparing and reclosing the NuDEL, flush the sheath (white hub) continuously with heparinised saline, the sheath is advanced to fully cover the mounted stent ensuring the delivering sheath is devoid of air.
2. Sterile preparation and draping of the access site should be performed according to standard practice or hospital protocol.
3. The entire system is advanced over the stiff guidewire into the desired location for implant. When advancing the NuDEL system over the wire, the tip of the wire must be controlled at all times.
4. After correct positioning of the stent, the sheath is pulled back to expose the stent; the large white hub assembly on the sheath should be pulled back to the double lines on the catheter shaft. The balloon / stent assembly may move distally a little during the process. Adjustment to achieve proper stent position is confirmed by a small injection of contrast through the sidearm of the sheath.
5. During inflation of the balloons, stability of the system versus the target lesion is determined by many factors such as guidewire stiffness, NuDEL stabilization, and blood flow. Consider using adequately stiff guidewires, and stabilize the outer sheath of the NuDEL to the groin.
6. Expand the stent initially by inflating the inner balloon using the inflation device until the inner balloon is fully expanded. The stent may be "repositioned" at this point by moving the NuDEL system. The unexpanded outer balloon and the expanded inner balloon hold the stent tightly against the NuDEL system. DO NOT deflate the inner balloon before expansion of the outer balloon. This could cause the stent to slip off the balloon catheter.
7. Confirming stent target position, inflate the outer balloon to the rated pressure. Do not exceed the manufacturer's balloon rated burst pressure.
8. Once the stent is expanded, deflate both balloons completely and rotate to ensure the stent is free and properly deployed. If there is a residual waist in the stent, expand only the outer balloon again, making sure not to exceed the rated burst pressure. A small contrast injection may be made through the sidearm of the sheath.
9. The deflated balloons are gently pulled back into the sheath, and the system removed. The result is confirmed by angiography.
10. Apply pressure to the insertion site according to standard practice or hospital protocol.

**NOTE: The diameter of the stent may be increased after placement by re-expansion with a larger diameter balloon. Do not exceed the maximum recommended expanded stent diameter (24mm). The stent will shorten during expansion (see Stent Foreshortening table).**

**RETURN OF EXPLANTED DEVICE:**

NuMED, Inc. is interested in obtaining recovered CP Stents. Place the explanted device in a container or vial immediately after excision. For further instructions on the return of an explanted device, contact RA Manager, NuMED, Inc. 2880 Main Street, Hopkinton, New York, 12965. Phone number: 315-328-4491.

**WARNING:**

NuMED stents are placed in the extremely hostile environment of the human body. Stents may fail to function for a variety of causes including, but not limited to, medical complications or failure of stent by fracture and embolization. In addition, despite the exercise of all due care in design, component selection, manufacture, and testing prior to sale, stents may be easily damaged before, during, or after insertion by improper handling, crimping or other intervening acts. Metal stents placed where there are extrinsic forces of compression, i.e. right ventricular outflow tract, are especially prone to fatigue fracture and embolization and should be avoided.

|                                 |
|---------------------------------|
| <b>Warranty and Limitations</b> |
|---------------------------------|

Stents and accessories are sold in an 'as is' condition. The entire risk as to the quality and performance of the stent is with the buyer. NuMED disclaims all warranties, expressed or implied, with respect to catheters and accessories, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. NuMED shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind NuMED to any representation or warranty with respect to catheters and accessories.

### CP Stent™ Foreshortening Chart

| Inflated Balloon Diameter | CP8Z16<br>(Stent length after expansion)<br>Percentage Shortening | CP8Z22<br>(Stent length after expansion)<br>Percentage Shortening | CP8Z28<br>(Stent length after expansion)<br>Percentage Shortening | CP8Z34<br>(Stent length after expansion)<br>Percentage Shortening | CP8Z39<br>(Stent length after expansion)<br>Percentage Shortening | CP8Z45<br>(Stent length after expansion)<br>Percentage Shortening |
|---------------------------|---|---|---|---|---|---|
| 12mm                      | (1.61) cm<br><b>2.8%</b>  | (2.18) cm<br><b>0.8%</b>  | (2.62) cm<br><b>4.4%</b>  | (3.23) cm<br><b>3.1%</b>  | (3.72) cm<br><b>1.9%</b>  | (4.17) cm<br><b>3.8%</b>  |
| 14mm                      | (1.54) cm<br><b>6.5%</b>  | (2.08) cm<br><b>5.4%</b>  | (2.56) cm<br><b>6.8%</b>  | (3.15) cm<br><b>5.4%</b>  | (3.66) cm<br><b>3.6%</b>  | (3.97) cm<br><b>8.4%</b>  |
| 15mm                      | (1.51) cm<br><b>8.5%</b>  | (2.02) cm<br><b>7.9%</b>  | (2.51) cm<br><b>8.6%</b>  | (3.10) cm<br><b>7.0%</b>  | (3.54) cm<br><b>6.6%</b>  | (3.94) cm<br><b>9.2%</b>  |
| 16mm                      | (1.48) cm<br><b>10.6%</b>   | (1.98) cm<br><b>10.1%</b>   | (2.45) cm<br><b>10.7%</b>   | (3.00) cm<br><b>9.8%</b>  | (3.48) cm<br><b>8.2%</b>  | (3.84) cm<br><b>11.4%</b>   |
| 18mm                      | (1.43) cm<br><b>13.7%</b>   | (1.89) cm<br><b>14.0%</b>   | (2.38) cm<br><b>13.3%</b>   | (2.88) cm<br><b>13.5%</b>   | (3.20) cm<br><b>15.6%</b>   | (3.71) cm<br><b>14.5%</b>   |
| 20mm                      | (1.32) cm<br><b>20.0%</b>   | (1.80) cm<br><b>17.9%</b>   | (2.30) cm<br><b>16.3%</b>   | (2.63) cm<br><b>20.9%</b>   | (2.96) cm<br><b>21.9%</b>   | (3.27) cm<br><b>24.7%</b>   |
| 22mm                      | (1.23) cm<br><b>25.4%</b>   | (1.67) cm<br><b>23.9%</b>   | (2.09) cm<br><b>24.0%</b>   | (2.46) cm<br><b>26.0%</b>   | (2.85) cm<br><b>25.0%</b>   | (3.15) cm<br><b>27.3%</b>   |
| 24mm                      | (1.05) cm<br><b>36.4%</b>   | (1.46) cm<br><b>33.8%</b>   | (1.91) cm<br><b>30.3%</b>   | (2.07) cm<br><b>37.9%</b>   | (2.27) cm<br><b>40.1%</b>   | (2.83) cm<br><b>34.9%</b>   |

### CP Stent™ Balloon Sizing Chart

| Inner Balloon Pressure (atm) | Stent ID (mm)              |                            |                            |                            |                            |                            |                            |
|------------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|
|                              | 12mm Diameter<br>RBP = 7.0 | 14mm Diameter<br>RBP = 6.0 | 16mm Diameter<br>RBP = 5.0 | 18mm Diameter<br>RBP = 4.0 | 20mm Diameter<br>RBP = 4.0 | 22mm Diameter<br>RBP = 3.0 | 24mm Diameter<br>RBP = 3.0 |
| 1                            | 4.40                       | 5.51                       | 6.57                       | 6.92                       | 7.71                       | 8.26                       | 8.92                       |
| 2                            | 4.85                       | 6.06                       | 7.03                       | 7.51                       | 8.45                       | 9.17                       | 9.84                       |
| 3                            | 5.18                       | 6.50                       | 7.40                       | 8.00                       | 8.93                       | 9.79                       | 10.51                      |
| 4                            | 5.56                       | 6.76                       | 7.78                       | 8.39                       | 9.44                       | 10.38                      | 11.23                      |
| 4.5                          |                            |                            |                            |                            |                            | <b>10.80</b>               | <b>11.76</b>               |
| 5                            | <b>5.89</b>                | <b>7.04</b>                | <b>8.09</b>                | <b>8.73</b>                | <b>9.98</b>                |                            |                            |
| Outer Balloon Pressure (atm) |                            |                            |                            |                            |                            |                            |                            |
| 1                            | 9.42                       | 10.84                      | 12.51                      | 13.92                      | 16.17                      | 18.44                      | 20.68                      |
| 2                            | 10.03                      | 11.61                      | 13.37                      | 15.14                      | 17.34                      | 19.87                      | 22.21                      |
| 3                            | 10.53                      | 12.19                      | 14.09                      | 16.36                      | 18.56                      | <b>21.81</b>               | <b>23.90</b>               |
| 4                            | 10.97                      | 12.76                      | 14.97                      | <b>17.69</b>               | <b>19.98</b>               |                            |                            |
| 5                            | 11.34                      | 13.44                      | <b>15.88</b>               |                            |                            |                            |                            |
| 6                            | 11.68                      | <b>13.89</b>               |                            |                            |                            |                            |                            |
| 7                            | <b>11.96</b>               |                            |                            |                            |                            |                            |                            |

The figures in bold face represent the stent ID @ Rated Burst Pressure.

**FOR ALL NUMED CATHETERS AN INFLATION DEVICE  
WITH PRESSURE GAUGE SHOULD BE USED.**

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G. van Wageningen B.V.