
# COVERED MOUNTED CP STENT™

## COARCTATION OF THE AORTA

Instructions for Use.....	2
Mode d'emploi .....	8
Istruzioni per l'uso.....	14
Gebrauchsanweisung.....	20
Instrucciones de utilización.....	27
Bruksanvisning .....	33
Brugsanvisning .....	39
Gebruiksaanwijzing .....	45
Instruções de utilização .....	51
Kullanım Talimatı .....	57

**CAUTION:** Federal (USA) Law restricts this device to sale by or on the order of a physician.

The ACC/AHA 2008 Guidelines for the Management of Adults With Congenital Heart disease recommends a yearly follow-up and additional imaging of the Coarctation site by CT or MRI every 5 years or less.

## Instructions for Use:

### INDICATIONS

Indicated for implantation in the native and/or recurrent coarctation of the aorta on patients with the following clinical conditions:

- Stenosis of the aorta resulting in significant anatomic narrowing as determined by angiography or non-invasive imaging, i.e. echocardiography, magnetic resonance imaging (MRI), CT scan;
- Stenosis of the aorta resulting in hemodynamic alterations, resulting in systolic pressure gradient, systemic hypertension or altered left ventricular function;
- Stenosis of the aorta where balloon angioplasty is ineffective or contraindicated;
- Stenosis diameter <20% of the adjacent vessel diameter. Stenosis that would present increased risk of vascular damage or disruption; or aneurysm associated with Coarctation of the aorta.

### DESCRIPTION

The Covered CP Stent is balloon expandable and intended for permanent implant. The Covered CP Stent is composed of heat treated 90% platinum/10% iridium wire that is arranged in laser welded rows with a "zig" pattern. The number of zigs in a row can be varied and will impact the strength of the stent as well as the eventual expanded diameter and percent stent shortening, while the number of rows will determine the unexpanded length of the stent. The Covered CP Stent has an ePTFE covering attached to the stent framework. This covering acts as a fluid barrier creating a fluid tight conduit through the stent length.

The NuMED BIB® (Balloon in Balloon) Catheter is a triaxial design catheter. Two lumens are used to inflate the balloons while one lumen is for tracking over a guidewire. The radiopaque platinum marker(s) are placed beneath the 'working area' of the balloon. The inner balloon is ½ of the outer balloon diameter and 1 cm shorter. Each balloon inflates to the stated diameter and length at a specific pressure. The balloon size is ±10% at rated burst pressure (RBP). The RBP is different for each size. Check the package label for the RBP. It is important that the balloon not be inflated beyond the RBP.

The purpose of the double balloon catheter is to apply an incremental inflation for the purpose of opening a vascular channel using a balloon expandable intravascular stent. The inner balloon provides initial expansion of the stent and also acts as a tool to hold the stent in place while the outer balloon is inflated. The outer balloon is then inflated securing the stent against the vessel wall.

### HOW SUPPLIED

Supplied sterilized by ethylene oxide gas. Sterile and non-pyrogenic if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

### CONTRAINDICATIONS

- Patients too small to allow safe delivery of the stent without compromise to the systemic artery used for delivery;
- Unfavorable aortic anatomy that does not dilate with high pressure balloon angioplasty;
- Occlusion or obstruction of systemic artery precluding delivery of the stent;
- Clinical or biological signs of infection;
- Active endocarditis;
- Known allergy to aspirin, other antiplatelet agents, or heparin;
- Pregnancy.

## **CP STENT WARNINGS**

- As with any type of implant, infection secondary to contamination of the stent may lead to aortitis, or abscess.
- The platinum/iridium stent may migrate from the site of implant.
- Over-stretching of the artery may result in rupture or aneurysm formation.
- The inflation diameter of the balloon used during stent delivery should approximate the diameter of the obstructive vessel and the intended implant site.

## **BIB STENT PLACEMENT WARNINGS**

- Do not exceed the RBP. An inflation device with pressure gauge is recommended to monitor pressure. Pressure in excess of the RBP can cause balloon rupture and potential inability to withdraw the catheter through the introducer sheath.
- Confirm that the distal end of the introducing sheath is at least 2.5cm back from the most proximal image markers before inflating the outer balloon. Failure to do so may stretch the outer tubing and severely hinder balloon deflation.
- Use two appropriate size inflation devices with pressure gauges for inflation.
- Do not advance the guidewire, balloon dilatation catheter, or any other component if resistance is met, without first determining the cause and taking remedial action.
- When the stent is crimped onto a balloon delivery catheter, the maximum balloon inflation pressure must not exceed the recommended inflation pressure specified in the manufacturer's instructions.
- This catheter is not recommended for pressure measurement or fluid injection.
- Do not remove the guidewire from the catheter at any time during the procedure except when the procedure has been completed.
- This device is intended for single use only. Do not resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of cross-contamination.

## **PRECAUTIONS**

- Use of an inflation device with pressure gauge is highly recommended during this procedure.
- Stents are delicate devices. Care should be exercised while handling to help prevent the possibility of breakage.
- The stent is rigid and may make negotiation through vessels difficult.
- Dilatation procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- Careful attention must be paid to the maintenance of tight catheter connections and by aspiration before proceeding to avoid air introduction into the system.
- The inflation diameter of the balloon used during stent delivery should approximate the diameter of the obstructive vessel and the intended implant site.
- Under no circumstances should any portion of the catheter system be advanced against resistance. The cause of the resistance should be identified with fluoroscopy and action taken to remedy the problem.
- If resistance is felt upon removal, then the balloon, guidewire and the sheath should be removed together as a unit, particularly if balloon rupture or leakage is known or suspected. This may be accomplished by firmly grasping the balloon catheter and sheath as a unit and withdrawing both together, using a gentle twisting motion combined with traction.
- Before removing the catheter from the sheath, it is very important that the balloon is completely deflated.
- Proper functioning of the catheter depends on its integrity. Care should be used when handling the catheter. Damage may result from kinking, stretching, or forceful wiping of the catheter.

## POTENTIAL COMPLICATIONS/ADVERSE EFFECTS

**NOTE:** Circumferential tear of the delivery balloon catheter prior to complete expansion of the stent may cause the balloon to become tethered to the stent, requiring surgical removal. In case of rupture of an adequately sized balloon after stent expansion, it can be withdrawn and a new balloon catheter exchanged over a guidewire to complete expansion of the stent. Cardiac catheterization carries certain risks. In addition, potential complications, and related adverse effects associated with implants include, but are not limited to:

- Stent migration
- Stent Fracture
- Vessel ruptures
- Minor hematoma
- Distal thromboemboli
- Death
- Intraluminal thrombus
- Hematoma requiring repair
- Bleeding requiring transfusion
- Cell necrosis at the site of implant
- Stent misplacement
- Pseudoaneurysm
- Stent misorientation
- Sepsis/infection
- AV fistula formation
- Transitory arrhythmia
- Cerebrovascular incident

## INSTRUCTIONS FOR USE

### Select Stent Size

1. Measure the length of the target stricture to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the stricture.
2. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent.  
Note: Should more than one stent be required, place the stent most distal from the puncture site first, followed by placement of the proximal stent in tandem.
3. Measure the diameter of the reference stricture and vessel proximal and distal to the target lesion to determine the appropriate size stent and delivery system.

### Preparation of Stent Delivery System

- Visually inspect the balloon/stent assembly to assure proper placement of the stent (MAY BE CONFIRMED BY FLUOROSCOPY).

### Stent Deployment

1. The assembly is advanced through the long delivery sheath and over the stiff guidewire into the desired location for implant. The supplied sleeve must be inserted in to the bleed back valve of the delivery sheath prior to insertion of the covered stent system. This is necessary to insure the covering remains intact. After the covered stent is completely beyond the valve, the protective sleeve can be pulled back to allow the valve to function properly.
2. After correct positioning of the stent, pull back on the sheath to expose the stent. Confirm proper stent position by a small injection of contrast through the sidearm of sheath or through a second catheter.
3. Expand the stent initially by inflating the inner balloon by twisting the locked pressure gauge inflation device until the inner balloon is fully expanded. One may "reposition" the stent at this point by moving the BIB® catheter. The unexpanded outer balloon and expanded inner balloon hold the stent tightly against the BIB® catheter.
4. Confirm positioning and inflate the outer balloon to rated diameter. Do not exceed the manufacturer's balloon rated burst pressure.

### Delivery System Withdrawal

1. Once the stent is expanded, deflate both balloons completely. If there is a residual waist in the stent, expand only the outer balloon again, making sure not to exceed RBP. Deflate the balloon and rotate the BIB® catheter to insure the stent is free and properly deployed.
  2. Remove the balloon catheter and confirm the result with angiography
- NOTE:** Diameter of the stent may be increased after placement by expanding with a larger diameter balloon. Do not exceed the maximum recommended expanded stent diameter.

**RETURN OF EXPLANTED DEVICE:**

NuMED, Inc. is interested in obtaining recovered CP Stents. Place the explanted device in a container or vial immediately after excision. For further instructions on the return of an explanted device, contact RA Manager, NuMED, Inc. 2880 Main Street, Hopkinton, New York, 12965. Phone number: 315-328-4491.

**WARNING:**

NuMED stents are placed in the extremely hostile environment of the human body. Stents may fail to function for a variety of causes including, but not limited to, medical complications or failure of stent by fracture and embolization. In addition, despite the exercise of all due care in design, component selection, manufacture, and testing prior to sale, stents may be easily damaged before, during, or after insertion by improper handling, crimping or other intervening acts. Metal stents placed where there are extrinsic forces of compression, i.e. right ventricular outflow tract, are especially prone to fatigue fracture and embolization and should be avoided.

<b>Warranty and Limitations</b>
---------------------------------

Stents and accessories are sold in an 'as is' condition. The entire risk as to the quality and performance of the stent is with the buyer. NuMED disclaims all warranties, expressed or implied, with respect to catheters and accessories, including but not limited to, any implied warranty of merchantability or fitness for a particular purpose. NuMED shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind NuMED to any representation or warranty with respect to catheters and accessories.

## CP Stent™ Foreshortening Chart

Inflated Balloon Diameter	CP8Z16 (Stent length after expansion) Percentage Shortening	CP8Z22 (Stent length after expansion) Percentage Shortening	CP8Z28 (Stent length after expansion) Percentage Shortening	CP8Z34 (Stent length after expansion) Percentage Shortening	CP8Z39 (Stent length after expansion) Percentage Shortening	CP8Z45 (Stent length after expansion) Percentage Shortening
12mm	(1.61) cm 2.8%	(2.18) cm 0.8%	(2.62) cm 4.4%	(3.23) cm 3.1%	(3.72) cm 1.9%	(4.17) cm 3.8%
14mm	(1.54) cm 6.5%	(2.08) cm 5.4%	(2.56) cm 6.8%	(3.15) cm 5.4%	(3.66) cm 3.6%	(3.97) cm 8.4%
15mm	(1.51) cm 8.5%	(2.02) cm 7.9%	(2.51) cm 8.6%	(3.10) cm 7.0%	(3.54) cm 6.6%	(3.94) cm 9.2%
16mm	(1.48) cm 10.6%	(1.98) cm 10.1%	(2.45) cm 10.7%	(3.00) cm 9.8%	(3.48) cm 8.2%	(3.84) cm 11.4%
18mm	(1.43) cm 13.7%	(1.89) cm 14.0%	(2.38) cm 13.3%	(2.88) cm 13.5%	(3.20) cm 15.6%	(3.71) cm 14.5%
20mm	(1.32) cm 20.0%	(1.80) cm 17.9%	(2.30) cm 16.3%	(2.63) cm 20.9%	(2.96) cm 21.9%	(3.27) cm 24.7%
22mm	(1.23) cm 25.4%	(1.67) cm 23.9%	(2.09) cm 24.0%	(2.46) cm 26.0%	(2.85) cm 25.0%	(3.15) cm 27.3%
24mm	(1.05) cm 36.4%	(1.46) cm 33.8%	(1.91) cm 30.3%	(2.07) cm 37.9%	(2.27) cm 40.1%	(2.83) cm 34.9%

## CP Stent™ Balloon Sizing Chart

Inner Balloon Pressure (atm)	Stent ID (mm)							
	12mm Diameter RBP = 7.0	14mm Diameter RBP = 6.0	15mm Diameter RBP = 5.0	16mm Diameter RBP = 5.0	18mm Diameter RBP = 4.0	20mm Diameter RBP = 4.0	22mm Diameter RBP = 3.0	24mm Diameter RBP = 3.0
1	2.75	3.22	3.49	3.75	3.94	4.02	4.20	4.28
2	2.85	3.32	3.59	3.85	4.36	4.13	4.33	4.50
3	5.85	6.91	6.89	7.79	8.54	9.20	10.16	10.57
4	6.12	7.00	7.02	7.95	8.71	9.63	10.40	11.08
4.5							<b>10.84</b>	<b>11.94</b>
5	<b>6.20</b>	<b>7.08</b>	<b>7.10</b>	<b>8.04</b>	<b>8.91</b>	<b>10.00</b>		
Outer Balloon Pressure (atm)								
1	10.73	13.08	13.45	14.87	16.85	17.91	20.52	22.79
2	10.86	13.27	14.16	15.10	17.06	18.38	21.46	23.95
3	11.15	13.50	14.55	15.68	17.64	19.42	<b>21.98</b>	<b>24.68</b>
4	11.33	13.68	14.88	15.93	<b>18.06</b>	<b>20.07</b>		
5	11.62	13.87	<b>15.06</b>	<b>16.19</b>				
6	11.80	<b>13.98</b>						
7	<b>12.04</b>							

\*This data is based on testing performed using the NuMED BIB® Stent Placement Catheter.

**The figures in bold face represent the stent ID @ Rated Burst Pressure.**

**FOR ALL NUMED CATHETERS AN INFLATION DEVICE  
WITH PRESSURE GAUGE SHOULD BE USED.**

<b>BIB DELIVERY CATHETER BALLOON DIAMETER AND INTRODUCER SIZE</b>	<b>REQUIRED INTRODUCER WITH BARE CP STENT</b>	<b>REQUIRED INTRODUCER WITH COVERED CP STENT</b>
<b>12MM (8F)</b>	<b>10F</b>	<b>12F</b>
<b>14MM (8F)</b>	<b>10F</b>	<b>12F</b>
<b>15MM (9F)</b>	<b>11F</b>	<b>12F</b>
<b>16MM (9F)</b>	<b>11F</b>	<b>12F</b>
<b>18MM (10F)</b>	<b>11F</b>	<b>14F</b>
<b>20MM (10F)</b>	<b>12F</b>	<b>14F</b>
<b>22MM (11F)</b>	<b>12F</b>	<b>14F</b>
<b>24MM (11F)</b>	<b>12F</b>	<b>14F</b>

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NuMED Canada Inc.